Eyelink Stimulus Presentation

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
官方编程软件：Experiment Builder

The image displays the SR Research Experiment Builder software interface, which is used for creating and running eye-tracking experiments. The interface is divided into several panels:

- Overview:** Shows a high-level flowchart of the experiment sequence.
- Structure:** A tree view showing the hierarchical organization of the experiment components, including DISPLAY_SCREEN, KEYBOARD, EL_BUTTON, TIMER, EL_CAMERA_SETUP, TRIAL, START, PREPARE_SEQUENCE, DRIFT_CORRECT, and RECORDING.
- Properties:** A table for configuring the selected component. The table has two columns: Property and Value.
- Diagram:** A visual representation of the experiment flow, showing the sequence of events and their timing.
- Navigator:** A panel at the bottom showing the current state of the experiment (BLOCK, TRIAL, RECORDING) and providing navigation controls.

The **Properties** table is as follows:

Property	Value
Label	RECORDING
Type	Sequence
Node Path	BLOCK.TRIAL.RECOR...
Time	
Record	☒
Recording Pause Time	20

The **Diagram** panel shows a flowchart with the following steps:

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

#1

Experiment Builder 2.3.38

SR Research Experiment Builder is a visual experiment creation tool that runs on Windows and macOS.

Recent / Future OS Releases

The release of Windows 11 and macOS BigSur and Monterey may involve changes that impact the functionality of our software. We are currently testing our software on these OS versions and will update their supported status as soon as our tests are complete.

Experiment Builder is backwards compatible (projects created in earlier versions can be opened and edited in later versions), and any existing deployed projects are standalone executables, so upgrading to the latest version will have no impact on data collection / ongoing research.

Please download the latest version from the links below.

Microsoft Windows 32-bit (x86) and 64-bit (x64) releases of Windows 7, 8, & 10

- Experiment Builder 2.3.38


macOS 10.12 Sierra to 10.15 Catalina

- Experiment Builder 2.3.38


Please note that changes in recent versions of macOS may show inconsistency in visual stimulus presentation timing. While macOS can be used to develop and debug experiments, if accurate display timing is critical we recommend deploying the task on a Windows PC for data collection.

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


Experiment Builder Graphical User Interface

Experiment Builder Life Cycle

- Experiment Design
- Building and Test-running Experiment
- Deploying Experiment
- Participant Data Set Randomization (optional)
- Data Collection
- Data Analysis

Working with Files——文件


- Creating a New Project
- Saving a Project
- Saving an Existing Project to a Different Directory
- Opening a Project
- Reopening a Recent Experiment Project
- Unlock/Lock a Project
- Packing an Experiment
- Unpacking an Experiment
- Experiment Node Exporting
- Experiment Node Importing
- Building an Experiment
- Cleaning an Experiment
- Test-running an Experiment from EB Application
- Deploying an Experiment
- Running an Experiment for Data Collection
- Converting Projects between Windows and macOS
- Running Experiments in EB 2.3

Designing an Experiment in Experiment Builder

- Hierarchical Organization of Experiments
- Experiment Graph: Flow Diagram
- Actions
- Triggers
- Other Node Types
- Sequence
- References and Equations
- Clearing Trigger Input Queue

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Dummy Mode

The screenshot displays the EYELINK software interface. The top panel, titled 'Structure', shows a tree view of devices with 'EYELINK' selected. Below it are tabs for 'Experiment', 'Components', and 'Devices'. The bottom panel, titled 'Properties', shows a table of configuration settings. The 'Dummy Mode' checkbox is currently unchecked. A yellow tooltip with the text 'false' is visible next to the 'Link Filter Level' property.

Property	Value
Tracker Address	100.1.1.1
Tracker Version	Eyelink 1000 Plus
Camera Mount	Current
Dummy Mode	☐
Data Processing	
Link Filter Level	NORMAL false
File Filter Level	HIGH
Eye Event Data	GAZE
Saccade Sensitivity	NORMAL
Eye Tracking	
Pupil Size	AREA
Fixation Update Interval	50
Fixation Update Accumulate	50
Auto Calibration Message	☒

Experiment Builder Project Check List

The screenshot displays the SR Research Experiment Builder software interface. The title bar reads "SR Research Experiment Builder (Document Version 2.3.1)". The interface is divided into a left sidebar with a tree view and a main content area on the right.

Left Sidebar (Tree View):

- Experiment Node Exporting
- Experiment Node Importing
- Building an Experiment
- Cleaning an Experiment
- Test-running an Experiment from EB Application
- Deploying an Experiment
- Running an Experiment for Data Collection
- Converting Projects between Windows and Mac OSX
- Running Experiments in EB 2.3
- Experiment Builder Graphical User Interface
 - Overview
 - Project Explorer Window
 - Graph Editor Window
- Application Menu Bar and Toolbar
- Designing an Experiment in Experiment Builder
 - Overview
 - Hierarchical Organization
 - Experiment Graph: Flow Diagram
 - Adding Components to an Experiment
 - Linking Experiment Components
 - Linking Rules
 - Actions
 - Triggers
 - Other Node Types
 - Sequence
 - References and Equations
 - Clearing Trigger Input Queue
 - Experiment Graph and Components
 - Overview
 - Graph Editing
 - Components
 - Screen Builder
 - Overview
 - Resources
 - Movement Patterns
 - Interest Areas
 - Resource Operations
- Data Source
- References
 - Overview
 - Using References
 - Entering in Values
 - Entering in References
 - Entering in Equations
 - Reference Manager
- EyeLink Data Viewer Integration
- Custom Class
- Examples
 - Experiment Builder Project Checklist
- Preference Settings
- Revision History

Main Content Area:

Experiment Builder Project Checklist (version 2.3.1)

Users may refer to the following checklist to prevent some common problems in creating and running an Experiment Builder project.

If this is an EyeLink experiment

- Is there a [recording sequence](#) in your EyeLink Enabled project?
When you check the "Record" property for a sequence in an EyeLink project, Gaze data will be recorded for the duration of that sequence (i.e., a command will be sent to the eye tracker to start recording when the sequence begins and a command will be sent to the eye tracker to stop recording when the sequence ends). This should typically only be enabled for the innermost sequence in the project. Given that recording sequences should only be used to control eye tracker recording and should not be used for looping purposes, the Iteration Count of recording sequences should always be set to 1 and there should never be a Data Source associated with a recording sequence.
- Does the project use the "hierarchical organization" concept of experiment design?
Employing a hierarchical design helps prevent redundancy and reduces overhead by allowing you to reuse the same trial structures to present multiple trials. For instance, sequences allow you to loop through the same sequence of nodes multiple times, while referencing values in the sequence's Data Source to change the properties of experimental events that need to change one trial to the next (e.g., to change the trials' images or sounds or to change the duration of timers from one trial to the next). Please see the Experiment Builder Video Tutorial Series for detailed instructions on how best to utilize hierarchical organization in your project. As a general rule, your project should have a single data source attached to the trial-level sequence. If your project has multiple data sources, or separate sequences/Data Sources for practice vs experimental trials, then it can make project editing and data analysis more difficult than they need to be.
- Has the 'Message' property of the critical nodes (e.g., [triggers](#) and [actions](#) used in the [recording sequence](#)) been filled out?
It is absolutely critical that you check through all the nodes in your recording sequence and make sure that all the important action and trigger nodes have appropriate text in the Message properties. Experiment Builder writes these messages to the EDF file to mark when these important events occur. Without these messages it may be impossible to align the gaze data with stimulus events accurately during data analysis or to use some Data Viewer features such as interest periods. It helps to use messages that will be readily interpretable at the analysis stage - e.g. STIMULUS_ONSET for the relevant display screen node, and TRIAL_END for the final display screen node etc. Please see the Experiment Builder examples for illustrations of the kinds of messages to send.
- Have you checked the Use for Host Display property of one (and only one) display screen node in your [recording sequence](#)?
You will not be able to see your stimuli on the Host PC recording unless you check this property for the most important display screen node in your recording sequence and set the "Draw To EyeLink Host" property of your prepare sequence to "image". Keep in mind only a single Display Screen action's contents can be transferred per trial. [Contents -> Experiment & Graph Components -> Components -> Actions -> Prepare Sequence]
- Has a [Recording Status Message](#) been written to the tracker screen to report progress of experiment testing?
The EyeLink Record Status Message property of a recording sequence allows you to pass a message to the Host PC to provide feedback to the experimenter during data collection. The message will appear in the bottom right of the Record screen on the Host PC, and can be used to provide information about the trial number, trial condition, etc. [Contents -> EyeLink to Data Viewer Integration -> Recording Status message]
- Has a [PREPARE_SEQUENCE](#) action been added before each [recording sequence](#)?
The prepare sequence node preloads stimulus resources (images, video, audio) for the upcoming recording sequence, transfers one display screen action's contents to the Host PC for experimenter feedback during the trial, and reinitializes the nodes of the upcoming sequence. It should be placed immediately before the drift correct action that precedes the recording sequence or immediately before the recording sequence if a drift correct is not used. Prepare sequences should not be placed at any higher levels in the experimental structure.
- Have you prepared [interest areas](#) for screen resources which you intended to analyze later?