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
SR Research论坛: <https://www.sr-support.com/>

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EyeLink® 1000 Plus User Manual

Desktop, LCD Arm, Tower, Primate
and Long Range Mounts

Remote, 2000 Hz and Fiber Optic
Camera Upgrades

Version 1.0.17

Copyright ©2013-2021, SR Research Ltd.

EyeLink is a registered trademark of SR Research Ltd.,

Oakville, Ontario, Canada

Table of Contents

1. Introduction	1
1.1 Supporting Documents	3
1.2 EyeLink 1000 Plus System Configuration	4
1.2.1 Host PC	4
1.2.2 Display PC	5
1.2.3 EyeLink 1000 Plus Camera Mount Configurations	6
1.3 System Specifications	9
1.3.1 Operational / Functional Specifications	9
1.3.2 Physical Specifications	10
2. EyeLink 1000 Plus Host Software	12
2.1 Web UI Interface	12
2.1.1 File Manager	12
2.1.2 Configuration Tool	17
2.1.3 Tracker Initialization Files	19
2.1.4 Running Web UI on a computer other than the host PC	20
2.2 Starting the Host Application	21
2.3 Modes of Operation	22
2.4 EyeLink 1000 Plus Host PC Navigation	23
2.4.1 Camera Setup Screen	24
2.4.2 Offline Screen	30
2.4.3 Set Options Screen	32
2.4.4 Calibrate Screen	40
2.4.5 Validate Screen	43
2.4.6 Drift Check/Drift Correct Screen	45
2.4.7 Output Screen	46
2.4.8 Record Screen	48
2.5 Status Panel	54
2.6 Mouse Simulation Mode	55
3. An EyeLink 1000 Plus Tutorial: Running an Experiment	56
3.1 The Camera Setup Screen	57
3.2 Participant Setup	58
3.2.1 Desktop Mount Participant Setup, Monocular	59

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

Reading & Language

Relevant Publications

Eye-tracking solutions for research in reading and spoken language comprehension.

Developmental

Relevant Publications

Understanding the development of cognitive, perceptual, social, and linguistic processes, from infants to the elderly.

fMRI & MEG

Relevant Publications

The world's most popular and effective solution for eye tracking in MRI/MEG environments.

EEG & fNIRS

Relevant Publications

Combine eye tracking with a range of other neurophysiological recording techniques for deeper insights.

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Clinical & Oculomotor

Relevant Publications


Research into oculomotor and cognitive dysfunction due to neuropsychiatric disorders.

Cognitive

Relevant Publications

Insights from eye tracking into attention, memory, decision making, and more.

Usability & Applied

Relevant Publications

Eye-tracking research out of the lab – solutions for usability and more.

Non-Human Primate


Relevant Publications

Fast, accurate, reliable, and non-invasive eye-tracking for non-human primates.

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All EyeLink Publications

All 11,000+ peer-reviewed EyeLink research publications up until 2022 (with some early 2023s) are listed below by year. You can search the publications library using keywords such as Visual Search, Smooth Pursuit, Parkinson's, etc. You can also search for individual author names. Eye-tracking studies grouped by research area can be found on the solutions pages. If we missed any EyeLink eye-tracking papers, please email us!

All years
▼

Search

11118 entries
▾ ▾ ▾
1 of 112
▹ ▹ ▹

2023

Chuanli Zang; Zhichao Zhang; Zhang Manman; Federica Degno; Simon P. Liversedge
Examining semantic parafoveal-on-foveal effects using a Stroop boundary paradigm Journal Article
In: Journal of Memory and Language, vol. 128, pp. 1–14, 2023.
[Abstract](#) | [Links](#) | [BibTeX](#)

Yuyang Zhang; Jing Yang; Zhisheng Edward Wen
Learners with low working memory capacity benefit more from the presence of an instructor's face in video lectures Journal Article
In: Journal of Intelligence, vol. 11, no. 5, pp. 1–14, 2023.
[Abstract](#) | [BibTeX](#)

Tania S. Zamuner; Theresa Rabideau; Margarethe McDonald; H. Henry Yeung
Developmental change in children's speech processing of auditory and visual cues: An eyetracking study Journal Article
In: Journal of Child Language, vol. 50, pp. 27–51, 2023.
[Abstract](#) | [Links](#) | [BibTeX](#)

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




Why thresholds are so important


Pupil and CR thresholds are greyscale values that can vary from 0 (black) to 255 (white).


You are basically telling the host how "black" or "white" something has to be to be considered as either pupil or CR.

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


Figure 2-4: Example Camera Setup Screen

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
Figure 1-1: Typical EyeLink 1000 Plus Configuration with a Desktop Mount

1.3 System Specifications


1.3.1 Operational / Functional Specifications

	Tower Mount /Primate Mount	Desktop and LCD Arm Mounts	
		Base System	Remote Tracking (Remote Camera Upgrade required)
Average Accuracy ¹	Down to 0.15° (0.25° to 0.5° typical)	0.25-0.5° typical	
Sampling rate ²	Monocular: 250,500,1000,2000 Hz Binocular: 250,500,1000,2000 Hz	Monocular: 250,500,1000 Hz Binocular: 250,500,1000 Hz	
End-to-End Sample Delay ³	1000 Hz: M = 1.97 ms, SD = 0.39 ms 2000 Hz: M = 1.34 ms, SD = 0.20 ms	500 Hz: M = 3.29 ms, SD = 0.58 ms 1000 Hz: M = 2.19 ms, SD = 0.30 ms	
Blink/Occlusion Recovery	1.0 ms @ 1000 Hz 0.5 ms @ 2000 Hz	2.0 ms @ 500 Hz 1.0 ms @ 1000 Hz	
Spatial Resolution ⁴	0.01°		
Noise with Participants ⁵	Filter (Off/Normal/High) 1000 Hz: 0.02°/ 0.01°/ 0.01° 2000 Hz (monoc): 0.03°/ 0.02°/ 0.01° 2000 Hz (binoc): 0.04°/ 0.02°/ 0.02°	Filter (Off/Normal/High) 500 Hz: 0.03°/0.02°/0.01° (25 mm lens) 0.06°/0.03°/0.01° (16 mm lens) 1000 Hz: 0.05°/0.03°/0.01° (25 mm lens) 0.08°/0.04°/0.02° (16 mm lens)	
Eye Tracking Principle ⁶	Dark Pupil - Corneal Reflection		
Pupil Detection Models	Centroid or Ellipse Fitting	Ellipse Fitting	
Pupil Size Resolution ⁵	0.1% of diameter	0.2% of diameter (16 mm lens) 0.1% of diameter (25 mm lens)	
Gaze Tracking Range	60° horizontally, 40° vertically	Customizable Default is 32° horizontally × 25° vertically	
Allowed Head Movements Without Accuracy Reduction	±25 mm horizontal or vertical	16 mm lens: 35×35 cm at 60 cm 40×40 cm at 70 cm 25 mm lens: 22×22 cm at 60 cm 25×25 cm at 70 cm	
Optimal Camera-Eye Distance	Tower: 48 cm Primate: 30 - 45 cm	40 - 70 cm	
Infrared Wavelength	Tower: 940 nm Primate: 910 /940 nm	850 to 940 nm	
Glasses Compatibility	Good	Excellent	Good
On-line Event Parsing	Fixation / Saccade / Blink / Fixation Update		
EDF File and Link Data Types	Gaze, Raw, and HREF eye position data/ Pupil size / Online events / Buttons / Messages / Digital inputs		
Real-Time Operator Feedback	Eye position gaze cursor superimposed on static image or position traces with camera images and tracking status.		

Specifications are preliminary and subject to change without notice.
¹ Measured with real eye fixations at multiple screen positions on a per subject basis.
² Availability depends on having the appropriate hardware and camera programming.
³ Time from physical event until first registered sample is available via Ethernet or Analog output. Optional data filter algorithm adds one sample delay for each filtering level.
⁴ Measured with an artificial eye.
⁵ Measured with real subject fixations.
⁶ Pupil-Only tracking mode is available for use in head fixed conditions.


1.3.2 Physical Specifications

Physical	Anodized aluminum enclosure.
GL (EyeLink 1000 Plus Camera)	Standard thread (¼"-20) centered on optical axis from 3 sides. M8 thread on front for DM and AM mounts. Power requirements: +12VDC, 800 mA for camera alone, 1.8A maximum when used with illuminator(s). Imaging rates up to 2000 fps.
OC (Fiber Optic Camera Head)	M3 tapped holes on body. Adapter for ¼"-20 mounting available. Supply: 3.6-5.7VDC, 700mA via LEMO connector and power harness. Imaging rates up to 2000 fps. Non-magnetic (< 0.5g of iron or nickel). Radio-quiet case and cable design for use in research environments.
FL-890, FL-940 (890/940 nm Fresnel Illuminator Module)	Standard thread (¼"-20) at center and sides. Adjustable focus and beamwidth. Wavelength: 890 or 940 nm. Supply: 3.6-5.7VDC, 1.1A via LEMO connector and power harness. Non-magnetic (< 0.5g of iron or nickel). Radio-quiet case and cable design for use in research environments.
PM-910, PM-940 (910/940 nm Illuminator Module)	Wavelength: 910 or 940 nm. Eye illumination level: less than 1 mW/cm ² at >200mm from illuminator. Illuminators powered from camera through supplied cables.
DM-850W, DM-850L, DM-890, DM-940, AM-890, AM-940 (DM and AM Series Illuminators)	Wavelength: 850 to 940 nm. Eye illumination level: less than 1 mW/cm ² at >450mm from illuminator. Powered from camera via integrated cables.
Camera Ethernet	Cabling: Unshielded CAT5e or CAT6 cables up to 30 meters in length. Requires host computer with supported Ethernet hardware capable of gigabit speeds.
Power Supply Specifications	GL camera: 12VDC, 2A external power supply with 2.5mm coaxial ("barrel") power connector (5.5 × 2.5 × 9.5mm). OC camera head and FL illuminator: 3.6-5.6VDC@2A minimum. When used with supplied power harness, DB-9 connector is required (contact SR Research for pinout). Power supply must have EN 60950, UL 950, CSA 22.2 No. 950, or other equivalent safety approval, with LPS or Class 2 rating.
Operating conditions	10°C to 30°C, 10%-80% humidity (non-condensing) For indoor use only.
Storage conditions	-10°C to 60°C, 10%-90% humidity (non-condensing). Allow to warm to room temperature before unpacking or use after storage at temperatures below 10°C to prevent condensation.
Safety Standards	IEC 60950-1:2005 (2nd Ed.) + Am 1:2009 IEC 60825-1:1993 + A1:1997 + A2:2001 IEC 62471:2006 (1st Ed.) UL 62368-1:2014 Ed.2 CSA C22.2#62368-1:2014 Ed.2 FDA 21 CFR laser products, under Laser Notice No. 50.
Fiber Optic Link (GL and OC modules)	Class 1 laser product, IEC 60825, CFR 21 850 nm, multimodal duplex fiber, LC connector.


2 Camera Connection

Using the **black** Ethernet cable, connect the Host PC Ethernet port labeled "EyeLink Camera" to the Ethernet port on the EyeLink 1000 Plus camera.

3 Host PC to Display PC Connection

Using the **blue** Ethernet cable, connect the Host PC Ethernet port labeled "Display PC Ethernet" to an Ethernet port on your Display PC. Configure that port with the following static IP address:

IP Address: 100.1.1.2

Subnet Mask: 255.255.255.0

Leave the default gateway and other settings blank.

See Chapter 9 of the EyeLink 1000 Plus Installation Guide for detailed instructions

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
Figure 9-2: Configure IP address on Windows 7

无法连接眼动仪: **Could not connect to tracker at 100.1.1.1**

“一万个主试会碰到一万次'**Could not connect to tracker at 100.1.1.1**', 但...”

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EyeLink 1000 Plus User Manual

EyeLink 1000 Plus Quick Start Guide

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

Figure 1-1: Suggested 1000 Plus System Layout

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Step A) Position the monitor so that it subtends no more than 32 degrees of visual angle horizontally and 25 degrees of visual angle vertically for the participant. The eye-to-monitor distance should be at least 1.75 times the display width to ensure that it falls within the trackable range.

Step B) Position the participant so that the eyes align with the top quarter of the monitor – adjust the chair and/or head support to ensure this alignment.

Step C) Position the eye tracker so that the distance from the top knob on the front of the Desktop Mount to the front of the chinrest is 50-55 cm. If using Remote Mode, make sure that the reported target distance on the Host PC is around 60 cm (55-60 cm is ideal for calibration).

Step D) Adjust the height of the eye tracker so that it is as high as possible without blocking the participant's view of the display.

Step E) Position the eye tracker so that its bottom knob is centered horizontally on the front of the monitor.

See Chapter 3 of the EyeLink 1000 Plus User Manual for detailed instructions

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7 Enter Screen Settings

Exit to the File Manager if the Host Application is running (press Ctrl-Alt-Q or click Offline and then Exit EyeLink).

From the File Manager, press the Configuration Button.

Configuration

From the Configuration screen press the Screen Settings button.

Follow the instructions that appear to enter the Screen Dimensions, Display Resolution, Eye-to-Screen Distance, and Camera-to-Screen Distance.

See Section 8.4 of the EyeLink 1000 Plus Installation Guide for detailed instructions


Figure 8-5: Updating Screen Dimensions

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